

GUERIDON PONT CHIRURGICAL

SURGICAL INSTRUMENT TABLE



αalphacoméd
MATÉRIEL MÉDICAL - VENTE & SAV

ÉQUIPEMENTS HOSPITALIERS

GUERIDON PONT CHIRURGICAL

SURGICAL INSTRUMENT TABLE

DESCRIPTION / DETAILS

La table pont FOURES est une déserte conçu pour faciliter l'ergonomie du poste de travail du chirurgien. Utilisé en bloc opératoire elle permet un passage intégral du plateau support au dessus de la table opératoire.

The FOURES surgical instrument table is a system designed to create the more ergonomics surgeon workstation. Used in operating theatre, it allow a free complete move over surgical table.

CARACTERISTIQUES / SPECIFICATIONS

- Hauteur Maximum 1,45 m – Minimum 1m
- Déplacement continu par manivelle
- Plateau acier inoxydable orientable
- Angle maximum 20 degrés avec butées
- Dimensions : 600 x 800 mm
- Pieds aluminium massif anodisé
- Quatre roulettes diamètre 100 mm, dont deux avec blocage
- Montants acier laqué
- Boîtiers d'angle aluminium anodisé

- Maximum height 1,45m – Minimum 1m
- Smooth telescopic move by crank
- Stainless steel tilting tray
- Maximum angle 20 degrees
- Dimensions: 600x800mm
- Solid anodised aluminium bases
- Four 100mm diameter castors two with brake
- Painted steel telescopic legs
- Anodised aluminium motion angle boxes

OPTIONS / OPTIONS

- Plateau bord relevé
- Plateau bord tombé
- Séparation centrale sur plateau

- Shelf with raised edge
- Shelf with fallen edge
- Shelf with central separation



ÉQUIPEMENTS HOSPITALIERS