

Products for all your Patient Positioning Needs

Robotic, Orthopedic, Gynecology, Urology,
Spine, Laparoscopy, Bariatric & advanced
surgical positioning accessories
to improve patient outcomes

NEW!

Provides optimal intraoperative positioning
for all open and arthroscopic limb procedures

See pages 12-13

Extra Surgical Hand | No Air or Power Needed | Attaches to Any OR Table Rail





The Future of Patient Positioning

the Better Choice

American Made—For Less

SchureMed is one of the world's leading manufacturers of operating room products specializing in surgical table accessories and positioners. All of our products are proudly made in the United States in our Abington, MA manufacturing plant and shipped worldwide.

Patient Positioning

Patient positioning is our priority and guides the design of all our products. We offer products for Gynecology, Urology, Orthopedic and General Surgery procedures. SchureMed offers a complete line of surgical table accessories for your positioning needs. Some of our featured products include our Great White Stirrups, Powered Shoulder Chair, and Schure Spine Frame.

Excellent Customer Service

At SchureMed our mission is to provide excellent customer service and to improve patient outcomes and safety for care providers. We are here to serve your needs with courtesy, respect and consideration. Our customer service representatives are knowledgeable, highly-trained professionals who will answer any questions you have about our products and help you begin using them with ease. Our sales representatives are ready to meet with you and are qualified to provide educational in-service sessions with your surgeons and staff upon request.

Expert Engineering

Our highly-trained engineers work closely with doctors, nurses and other members of the medical community to incorporate user feedback into our products. The clinical feedback we receive allows us to produce high quality positioning devices that fit the needs of our customers.

High Quality Products at a Reasonable Price

Our in-house design and manufacturing team allow us to produce superior products at a more reasonable price than our competitors. Please visit our website at www.schuremed.com.



TABLE OF CONTENTS



Armboards	Patient Transfer.....	43
Armboard, Adjustable Height.....	Shoulder	
Armboard, Extra-Wide.....	Complete Shoulder Systems	15, 17
Armboard, Holder.....	Disposable Face Masks.....	15, 17, 18
Armboard, Incline Pad.....	Disposable Hand/Forearm Wraps	15, 17, 18
Armboard, Radiolucent.....	Disposable Head & Chin Straps.....	15, 17, 18
Armboard, Standard.....	Dolly, Shoulder Chair	15, 17
Floor Mat	Dolly, Shoulder Assistant	18
Gynecology/Laparoscopy/Urology	E-Z Lift Beach Chair	16-17
Candy Cane Stirrups	Lateral Shoulder Holder.....	19
Disposable Clear View Leg Drapes	Powered Beach Chair	14-15
Disposable Uro-Trapper Drain Bags.....	Schure Loc™ XPS.....	12-13
Dolly, Great White Stirrup.....	Shoulder Assistant	18
Great White Maxima Stirrups	Spine	
Great White Platinum Stirrups.....	Dolly, Schure Spine Frame	27
Great White Premium Stirrups	Schure Spine Frame.....	27
Great White Robotic Stirrups	Table Accessories	
Knee Crutches	Anesthesia Screen	35
Lithotomy Stirrups	Bariatric Nissen Strap.....	33
Pediatric Stirrups.....	Bariatric Restraint Straps.....	33
Split Leg Foot Boards.....	Disposable A-Line Positioners.....	39
Split Leg Stirrups	Extremity Positioner	36
Stirrup Wall Rack	Foot Extension	32
Uro-Trapper System.....	Irrigation Towers.....	38
Hand/Wrist/Elbow/Arm	Lateral Braces.....	34
Elbow Arthroscopy Support.....	Lateral Positioners	34
Elbow Surgical Positioner.....	Leg Prepper	35
End Rest Procedure Table	Narrow Headrest.....	32
Hourglass Procedure Table	Nissen Straps	33
Major Procedure Table	Restraint Strap	33
Minor Procedure Table.....	Safety Draw	38
Multi-Axis Arm Positioner.....	Saphenous Vein Harvesting Pad.....	36
Wrist Tower	Schure Strap VELCRO® brand Rolls.....	39
Knee	Shoulder Supports	33
Arthroscopic Leg Holder.....	Table Width Extenders.....	32
Disposable TKR Boot Liners.....	Toboggan Arm/Leg Guards	35
Schure Foot for Total Knee Replacement.....	Universal Prepper.....	37
Schure Grip Leg Holder	Table Clamps	
TKR Positioner for Total Knee Replacement... 28-29	CamLoc Clamp.....	40
Well Leg Holder	Deluxe Rail Clamp.....	41
Pads	Schure Socket	41
Armboard Pads.....	Schure Socket XL.....	40
Coated Positioning Pads	Simple Clamp	41
Disposable Positioning Pads	SpringLoc Clamp	40
Gel Pads	Universal Rail Clamp	41
Stretcher Pads		
Surgical Table Pads		

All products are not made with natural rubber latex.

600 lb. (272 kg)
Patient Wt. Capacity

Great White Platinum Stirrups

Exclusive
axial rotation
from the handle
for superior
patient positioning



Exclusive
3-axis control

Patent Pending

GYNECOLOGY/LAPAROSCOPY/UROLOGY

Great White Platinum Stirrups

- Exclusive 3-axis control from the handle
- Surgical weight capacity up to 600 lbs
- Boots free-float throughout lithotomy range
- Clam shell stirrup boot liners contain the foot, ankle and calf for maximum patient security and comfort
- Lithotomy and length indicators assist with accurate positioning of legs

The Great White Platinum Stirrups provide safe and quick positioning assistance for improved surgical site access. This intra-operative, adjustable medical stirrup for gynecological, laparoscopic, and urological applications offers exceptionally precise placement, effectively minimizing the possibility of nerve damage. Exclusive axial rotation, high/low lithotomy and abduction/adduction from the handle improves patient positioning.

Simply squeeze the trigger for smooth and accurate patient adjustment without jeopardizing sterility or patient comfort. Once the trigger is released, the stirrup is firmly locked for optimum safety and surgical performance. The mounting blade design quickly and easily inserts into side-rail clamps. The raised lateral boot fin system protects the patient from superficial nerve damage.

The enhanced range of motion allows for -35° for low lithotomy, to 84° high lithotomy; abduction -9° to abduction 25° and also externally rotates out 20°. The Great White Platinum Stirrups are delivered fully assembled with included boot liners. Two required SpringLoc Clamps sold separately.

#800-0007 Great White Platinum Stirrups

Required Accessory: Purchase Two

#800-0127 SpringLoc Clamp

Accessories

#800-0075 Dolly, Great White Stirrup

#800-0284 Great White Stirrup Wall Rack

#508-0113 Reusable Boot Pads

#508-0350 Yellow Fin Replacement Pad

Organize your Stirrups & Accessories



#800-0075 Great White Stirrup Dolly



#800-0284 Great White Stirrup Wall Rack



Great White® Premium Stirrups

- Exclusive axial rotation from the handle for superior patient positioning
- Lift assist allows easy control with minimal effort
- Flexible boot accommodates larger patients
- 400 lb patient weight capacity
- Better patient ergonomics
- Rotates on three axes
- Boot pads are included

High quality, feature-rich, and priced right, the Great White Premium Stirrups are exactly what your surgical team requires. The enhanced range of motion allows for -35° for low lithotomy, to 84° high lithotomy; abduction -9° to abduction 25° and also externally rotates out 20°.

#800-0301 Great White® Premium Stirrups, 11"W x 12"D x 36"H, 9.5 lbs

Required Accessory

#800-0127 SpringLoc Clamp (2 required)

Optional Accessories

- #800-0160** Disposable Clear View Leg Drapes, Sterile, 10/cs
- #800-0075** Dolly, Great White® Stirrup
- #800-0284** Great White Stirrup Wall Rack
- #508-1354** Replacement Premium Boot Pad (w/o Fin)
- #508-0113** Replacement Premium Boot Pad (w/Fin)

Allow visibility of patient's leg during intra-operatively adjustments with Disposable Clear View Drapes (left). Organize your accessories with the Stirrup Dolly (center) and Great White Stirrup Wall Rack (right).



Exclusive axial rotation from the handle for superior patient positioning

400 lb. (272 kg) Patient Wt. Capacity

Exclusive 3-axis control

Patent Pending



Split Leg Stirrups

- Superior split leg stirrups for lower extremity procedures
- Use for Vascular & Orthopedic procedures and mid-line access to abdomen

The Split Leg Positioner provides a safe system for positioning patients for laparoscopic procedures as well as lower extremity orthopedic and vascular procedures. The Split Leg Positioners are radiolucent and lightweight. They are designed with excellent access to the operative site. The mounting blade design quickly and easily inserts into side-rail clamps.

#800-0318 Split Leg Stirrups

Required Accessory: Order Two

#800-0336 Split Leg Clamp

Optional Accessories

- #800-0075** Dolly, Great White® Stirrup
- #800-0284** Great White Stirrup Wall Rack
- #800-0324** Split Leg Foot Boards



Great White Robotic Stirrups

- 600 lb patient weight capacity
- Rotates on three axes
- Flexible boot accommodates larger patients
- Lift assist allow easy control with minimal effort

The mounting blade design quickly and easily inserts into side-rail clamps. The Great White Robotic Stirrups' range of motion is 55° low lithotomy to 60° high lithotomy. Adduction -9° and abduction +25° and also externally rotates out 20°, with the control to lock the medical stirrup for positioning anywhere in between. Two required SpringLoc Clamps sold separately.

#800-0330 Great White Robotic Stirrups w/Fin
#800-0329 Great White Robotic Stirrups w/o Fin

- #508-0113** Replacement Boot Pad (w/Fin)
- #508-1354** Replacement Boot Pad (w/o Fin)

Required Accessory: Order Two
#800-0127 SpringLoc Clamp



Optional Accessories

- #800-0075** Dolly, Great White® Stirrup
- #800-0284** Great White Stirrup Wall Rack



Pediatric E-Z Lift Stirrups

The Pediatric E-Z Lift Stirrups are as easy-to-use as its forerunner. Exclusive axial rotation from the handle for superior patient positioning. Boot adjusts with a single clamp which simplifies patient positioning. Boot pads included.

- #800-0231 PUPS, 160 lb (73 kg) patient capacity
- #800-0232 KIDS, 250 lb (113 kg) patient capacity

#800-0238 Pediatric E-Z Lift Stirrups Bundle
Includes 1 pair of each stirrup

Required Accessory: Order Two
#800-0127 SpringLoc Clamp

Optional Accessories

- #800-0259 Dolly, Pediatric Stirrup
- #508-1295 PUPS Reusable Boot Pads
- #508-1296 KIDS Reusable Boot Pads



Exclusive axial rotation from the handle for superior patient positioning

- E-Z Lift feature easily control the patient's leg when moving it into the operative position
- Floating boot automatically adjusts when repositioning the pediatric boot stirrup
- Range: Abduction +25°, Adduction -9°, Lithotomy range of +84° to -35° & externally rotates out 20°

Pediatric Lithotomy Stirrups

The Pediatric Lithotomy Stirrups are a comfortable boot stirrup. Specifically designed with the pediatric patient in mind. Pediatric stirrups attach to side rail of any OR table using the Schure Socket XL (sold separately). Boot pads included.

- #800-0233 PUPS, 160 lb (73 kg) patient capacity
- #800-0234 KIDS, 250 lb (113 kg) patient capacity

#800-0237 Pediatric Lithotomy Stirrups Bundle
Includes 1 pair of each stirrup

Required Accessory: Order Two
Schure Socket XL

- #800-0134 US #800-0200 UK
- #800-0244 EU #800-0201 DEN

Optional Accessories

- #508-1295 PUPS Reusable Boot Pads
- #508-1296 KIDS Reusable Boot Pads



- Pediatric boot stirrup design encapsulates patient's lower leg
- Floating boot automatically adjusts when repositioning the pediatric boot stirrup

800 lb. (363 kg)
Patient Wt. Capacity



Exclusive 3-axis control

Patent Pending



Fold-over design



Flexible boot

Great White Maxima Stirrups

- Exclusive axial rotation from the handle for superior patient positioning
- A wide boot and a stronger gas cylinder make this stirrup perfect for your heaviest patients
- 800 lb (363 kg) surgical patient weight capacity
- Boot free-floats throughout each range reducing the risk of superficial nerve injury
- Lithotomy and length indicators help to position and reposition the legs
- Extended lateral boot design helps reduce the risk of peroneal nerve injury

Increase your surgeons' operating room performance with the Great White Maxima Stirrups which permit simple adjustment of lithotomy and abduction intraoperatively. Exclusive axial rotation from the handle for superior patient positioning. Great White Maxima Stirrups range of motion is -35° low lithotomy to 84° high lithotomy. Adduction -9°, abduction 25° and also externally rotates out 20°, with the control to lock the stirrup anywhere in between. The Great White ships fully assembled. Boot pads are included. Two required SpringLoc Clamps sold separately.

#800-0245 Great White Maxima Stirrups

Required Accessory: Purchase Two
#800-0127 SpringLoc Clamp

Optional Accessories

- #800-0075 Dolly, Great White Stirrup
- #800-0284 Great White Stirrup Wall Rack
- #508-0559 Replacement Bariatric Boot Pads



Uro-Trapper System

- Flexible drape support provides unobstructed surgeon access to surgical site and springs back to original orientation to allow for optimal fluid drainage
- Mounts effortlessly to any surgical table with integral clamps

A flexible system that attaches to any general surgical table to contain fluids for hysteroscopy or cystoscopy procedures.

#800-0096 Uro-Trapper System



Lithotomy Stirrups

- Use for long or short procedures
- Soft boot pads encapsulate feet
- 350 lb (158 kg) patient weight capacity
- Features extended lateral boot support to prevent superficial nerve damage and unwanted leg rotation
- Schure Socket XL sold separately

#800-0049 Lithotomy Stirrups

#800-0134 Schure Socket XL



Uro-Trapper Drain Bags

- Flexible hose is compatible with all suction systems for accurate fluid measurement
- Durable, drain bags are not made with natural rubber latex

Open mesh filter allows for rapid drainage of fluids. Tube is 64"L x 7/8"D (163cm x 2cm) and 20" (51cm) compressed.

Disposable Uro-Trapper Drain Bags

#800-0098 Sterile, 10/cs

#800-0095 Non-Sterile, 10/cs

Knee Crutches

- Simple to set up and position the patient
- Knee Crutches adjust vertically from 12"-19" (31cm-48cm) and rotate 360° to provide precise positioning with a simple ergonomic locking handle
- Constructed of the highest quality stainless steel
- Extremely soft and comfortable crutch liners prevent superficial nerve damage to the patient. Mount crutches simply with Schure Socket XL (sold separately)

#800-0041 Knee Crutches

#800-0134 Schure Socket XL



Disposable Clear View Leg Drapes

The Disposable Clear View Leg Drape allows visibility of patient's leg during intraoperative adjustments. A blue modesty panel is now incorporated within the drape. Packaged sterile and not made with natural rubber latex. Provides continuous monitoring of patient's legs without breaking sterile field.

#800-0160 Disposable Clear View Leg Drapes, Sterile, 10/cs

Candy Cane Stirrups

- Adjust vertically from 28"-44" (71cm-112cm) with an 8" (20cm) leg clearance
- 18"L x 1 1/2"W (46cm x 4cm) double-looped, cotton ankle strap
- Lock stirrup height by simply twisting the ergonomic handle
- Connect stirrups to OR table side rail with Schure Socket XL (sold separately)

#800-0012 Candy Cane Stirrups


#800-0134 Schure Socket XL

#790-0006 Candy Cane Strap Replacement

#512-0003 Replacement Hooks



Schure Loc™ XPS



Easy as 1, 2, 3

- Squeeze and hold lever to release all device joints
- Position limb in exact position
- Release lever to securely lock into position

The Schure Loc™ XPS combines strength and flexibility to provide optimal intraoperative positioning for all open and arthroscopic limb procedures.

State-of-the-art patient positioning device improves OR staff performance

Patent Pending

Schure Loc™ XPS

- Securely and safely holds patient arm in any desired position during arthroscopic surgery
- Compact and lightweight
- No foot pedals, hoses, hydraulics or power cords
- Ships assembled and ready for immediate use
- Single lever releases all joints for positioning and relocks when released
- Decreases surgeon fatigue and improves intraoperative control of surgical limb
- Schure Loc™ XPS reduces staffing needs providing a more cost-effective procedure
- Schure Loc™ XPS Shoulder Kit conveniently includes a sterile drape, foam hand/forearm holder and self-adherent wrap
- 500 lb weight capacity

Schure Loc™ XPS currently outperforms similar equipment for shoulder surgery in both the beach chair and lateral position. The Schure Loc™ XPS is the perfect companion piece for any beach chair positioner—easily attaching to any OR table rail while freeing your surgical assistant to do much more than steady the patient's limb. Get the most from your team by adding this high-quality, easy-to-use device to your shoulder surgery equipment.

#800-0280 Schure Loc™ XPS

Optional Accessories

#800-0281 Disposable Schure Loc™ XPS Sterile Kit, 10/cs

Includes Sterile Forearm Wrap, Sterile Self-Adherent Coban™ and Sterile Drape (fits whole arm).

#790-0039 Schure Loc™ Adapter

#790-0040 Schure Loc™ Forearm Attachment



Provides excellent positioning and surgical access for beach chair, lateral shoulder and supine elbow procedures.

Completely eliminates manual positioning with push-button ease

Powered Beach Chair

NEW! Double Ball Joint Head System allows complete control of the patient's head with one knob. Multiple articulating joints allow the head positioner to safely and efficiently control the most challenging patient anatomies, even super kyphosis patients.

Lateral counter-traction brace adjusts to secure patient positioning.



Head system slides vertically and then locks tight.



Patented

Powered Beach Chair

- 500 lb (227 kg) patient weight capacity
- Motion-control technology eliminates manual patient positioning
- Lateral Brace attaches to chair and not table side rail
- Pull-away panels providing excellent surgical exposure
- Double-ball joint head rest tilts 30° in every direction
- Five-piece Deluxe pad set included

The SchureMed Powered Beach Chair is the only chair on the market with a self-contained motor that does not have to rely on the surgical table for power. The double-ball joint head system safely secures patient's head; with even your kyphosis patients. Chair adjusts from 0-90° in 8 seconds in a slow controlled movement. 2" (5cm) Deluxe pad set provides patient comfort. Motion-control technology eliminates manual patient positioning with a simple push button hand control and supports patients up to 500 lbs (227 kgs). Attaches to ANY operating room table with our integrated adjustable mounting clamps. Lateral Brace included.

#800-0004 Powered Beach Chair

Required Accessories: Select One Model

#800-0164 Disposable Face Masks, 12/cs

#800-0039 Disposable Head & Chin Straps, 12/cs

Optional Accessories

#800-0074 Dolly, Shoulder Chair

#800-0050 Multi-Axis Arm Positioner

#800-0134 Schure Socket XL



Disposable Face Mask



Disposable Head & Chin Strap



Multi-Axis Arm Positioner



Schure Socket XL



39"H x 27"W x 12"D

Shoulder Chair Dolly



Complete Shoulder System

Includes the Schure Loc™ XPS, Disposable Head & Chin Straps, Multi-Axis Arm Positioner, Schure Socket XL, Disposable Schure Loc™ XPS Sterile Kit and Shoulder Chair Dolly.

Choose Powered Beach Chair or E-Z Lift Beach Chair

#800-0293 Complete Shoulder System w/Schure Loc™ XPS & Powered Beach Chair

#800-0294 Complete Shoulder System w/Schure Loc™ XPS & E-Z Lift Chair

E-Z Lift Beach Chair

NEW! Double Ball Joint Head System allows complete control of the patient's head with one knob. Using multiple articulating joints allow the head positioner to safely and effectively control the most challenging patient anatomies, even super kyphosis patients.



E-Z Lift Beach Chair

- Reduces effort when lifting patients compared to manual systems
- Double-ball joint head system safely secures patient's head

The E-Z Lift Beach Chair provides a lightweight and cost-effective option to our fully powered shoulder chair. This system employs a powerful piston to lift the patient. Adjusts to fit any OR table with our integrated adjustable mounting clamps. 500 lb (227 kg) patient weight capacity. Articulation range is 0°–90°. Lateral Brace included.

#800-0142 E-Z Lift Beach Chair

Required Accessories: Select One Model

#800-0164 Disposable Face Masks, 12/cs

#800-0039 Disposable Head & Chin Straps,

Optional Accessories

#800-0074 Dolly, Shoulder Chair

#800-0050 Multi-Axis Arm Positioner

#800-0134 Schure Socket XL



Disposable Face Mask



Disposable Head & Chin Strap



Shoulder Chair Dolly



Schure Socket XL



Multi-Axis Arm Positioner



Complete Shoulder System

Includes the Shoulder Assistant, Disposable Head & Chin Straps, Multi-Axis Arm Positioner, Shoulder Assistant Dolly, Schure Socket XL and the Disposable Hand/Forearm Wraps and Shoulder Chair Dolly

Choose E-Z Lift Beach Chair or Powered Beach Chair

#800-0130 Complete Shoulder Assistant Powered Beach Chair Shoulder System

#800-0295 Complete Shoulder Assistant E-Z Lift Chair Shoulder System



Shoulder Assistant

- Save money not having to buy expensive drapes
- Surgeon foot control provides quick and precise intraoperative positioning
- Easy-to-assemble
- Shoulder Assistant Dolly sold separately
- Unlimited surgical site access and positioning

The SchureMed Shoulder Assistant combines strength and adaptability to provide optimal intraoperative positioning for shoulder surgery in the beach chair position. Staff efficiency is improved as the assistants are freed from holding the limb during surgery. Designed to maintain a low profile while delivering exceptional control and positioning capabilities. Mounts onto the side of any standard operating table and provides a non-sterile base for which sterile components are inserted.

#800-0088 Shoulder Assistant

Required Accessory
Disposable Sterile Hand/Forearm Wraps securely holds the hand and forearm during the operative procedure. Quick release for range of motion.



#800-0113 Disposable Sterile Hand/Forearm Wraps, 12/cs

Optional Accessories

#800-0039
Disposable Head & Chin Straps, 12/cs

#800-0164
Disposable Face Masks,

#800-0140
Dolly, Shoulder Assistant



#800-0039



#800-0164



#800-0140



Eliminates use of sandbags

Lateral Shoulder Holder

This weightless shoulder system allows for traction settings from 0-25 lbs (0-11 kgs) by simply turning the traction adjustment knob. Internal and external rotation can be easily controlled by the surgeon's preference, and locked in place. It simply attaches to the surgical table side rail. No need for sand bags!

Required CamLoc Clamp sold separately.

#800-0169 Lateral Shoulder Holder

Required Accessory

CamLoc Clamp: Order One

#800-0271 US

#800-0273 UK

#800-0272 EU

#800-0274 DEN

#800-0113 Disposable Hand/Forearm Wraps, Sterile, 12/cs

Wrist Tower

- *Simply dial in your tension...eliminates use of sand bags*
- *All components are autoclavable*
- *Easily attaches to side rail of surgical table with CamLoc Clamp (sold separately)*
- *Improves productivity... eliminates the need for assistant to hold extremity*

The Wrist Traction Configuration offers a variety of positioning capabilities. You can apply traction by simply dialing the knob from 0-25 lbs (0-11 kgs). Pronation and Supination can be easily controlled by the surgeon's preference and locked in place.

#800-0135 Wrist Tower

Required Accessory: Order One CamLoc Clamp

#800-0271 US **#800-0273** UK
#800-0272 EU **#800-0274** DEN

Optional Accessory

#512-0015 Finger Traps, Limited Reuse
Please call customer service for product pricing and specifications.



Attach to table with CamLoc Clamp

Elbow Surgical Positioner

- *All components are autoclavable*
- *Improved productivity...no need for assistant to hold extremity*

The Elbow Surgical Positioner combines ease of use and adaptability to provide intraoperative positioning for supine elbow surgery. Quick release allows for intraoperative check for range of motion.

#800-0275 Elbow Surgical Positioner

Required Accessory

CamLoc Clamp: Order One

#800-0271 US **#800-0273** UK
#800-0272 EU **#800-0274** DEN

#800-0113 Disposable Sterile Hand/Forearm Wraps, 12/cs



End Rest Procedure Table

Choose Carbon Fiber or Phenolic

- Easy-to-use durable surgical table
- Carbon fiber technology produces incredibly solid components that are completely invisible to x-rays
- Adjusts up or down to be level with the OR table
- End rail accommodates traction attachment and locks to side rail in seconds
- Adjustable leg height with twist handle locking system

Table measures 16"W (41cm) at narrow end and 34"L x 26" (86cm x 66cm) at wider end. Deluxe foam pad included.



Hourglass Procedure Table

Choose Carbon Fiber or Phenolic

- Easy-to-use durable surgical table
- Carbon fiber technology is lighter, invisible to x-rays and provides a clearer picture during x-ray
- Adjusts up or down to be level with the OR table
- End rail accommodates traction attachment and locks to side rail in seconds
- Adjustable leg height with twist handle locking system

Table measures 15"W (38cm) at middle and 32"L x 23"W (84cm x 58cm) at table ends. Deluxe foam pad included.



Ordering Information

End Rest Major Procedure Table

- #800-0029 End Rest Major Procedure Table, Carbon Fiber
- #800-0035 End Rest Major Procedure Table, Phenolic
- #508-0101 Replacement Deluxe Foam Pad

Hourglass Major Procedure Table

- #800-0031 Hourglass Major Procedure Table, Carbon Fiber
- #800-0040 Hourglass Major Procedure Table, Phenolic
- #508-0106 Replacement Deluxe Foam Pad

Two Clamps Required for Either Product: Order Two

- Universal Clamp**
- #800-0085 US
- Deluxe Rail Clamp**
- #800-0205 EU
- #800-0206 UK
- #800-0207 DEN



Major Procedure Table

- Choose Carbon Fiber or Phenolic*
- Easy-to-use durable surgical table
 - Carbon fiber technology produces incredibly solid components that are completely invisible to x-rays
 - Adjusts up or down to be level with the OR table
 - End rail accommodates traction attachment and locks to side rail in seconds
 - Adjustable leg height with twist handle locking system

Table measures 15"W x 34"L (38cm x 86cm).
1" Deluxe foam pad & Leg included.

Standard Armboard

- Simply place armboard on OR table side rail and it snaps locked tight
- Armboard offers a full range of motion, swings 180°
- Simply pull finger trigger at end of armboard and move to required position. Release trigger and it locks tight. To remove, just pull finger trigger and lift off side rail.


- Armboard, Standard**
- | | | | |
|-----------|-----------------------|-----------|-----|
| #800-0008 | US | #800-0189 | UK |
| #800-0188 | EU | #800-0190 | DEN |
| #508-0071 | 2" (5cm) Deluxe Pad | | |
| #508-0151 | 2" (5cm) Softcare Pad | | |



Pad sold separately



New trigger for quick & easy removal



Minor Procedure Table

- Choose Carbon Fiber or Phenolic*
- Easy-to-use durable surgical table
 - Carbon fiber technology produces incredibly solid components that are completely invisible to x-rays
 - Adjusts up or down to be level with the OR table
 - Perfect for all arm and hand procedures
 - No leg attachment

Table measures 15"W x 30"L (38cm x 76cm).
1" Deluxe foam pad included.

Adjustable Height Armboard

- This armboard adjusts vertically to be level with OR table, either with or without X-ray tops
- To adjust armboard vertically, simply place armboard on side rail, loosen handle and raise to correct table height and tighten handle
- Eliminates multiple size pads

- Armboard, Adjustable**
- | | | | |
|-----------|-----------------------|-----------|-----|
| #800-0009 | US | #800-0192 | UK |
| #800-0191 | EU | #800-0193 | DEN |
| #508-0071 | 2" (5cm) Deluxe Pad | | |
| #508-0151 | 2" (5cm) Softcare Pad | | |



Pad sold separately

Ordering Information

Rectangle Major Procedure Table

- #800-0024 Carbon Fiber
- #800-0056 Phenolic
- #508-0073 Replacement 1" Deluxe Foam Pad

Rectangle Minor Procedure Table

- #800-0027 Carbon Fiber
- #800-0026 Phenolic
- #508-0072 Replacement 1" Deluxe Foam Pad

Two Clamps Required for Either Product: Order Two

- Universal Clamp**
- #800-0085 US
- Deluxe Rail Clamp**
- #800-0205 EU
 - #800-0206 UK
 - #800-0207 DEN

Extra-Wide Armboard

- Armboard, Extra-Wide**
- 11"W x 26½"L (28cm x 67cm) to offer a bigger surface area when necessary

- | | | | |
|-----------|----|-----------|-----|
| #800-0018 | US | #800-0195 | UK |
| #800-0194 | EU | #800-0196 | DEN |

- Armboard, Extra-Wide, Adjustable**
- This armboard adjusts vertically to be level with OR table, either with or without X-ray tops

- | | | | |
|-----------|----|-----------|-----|
| #800-0017 | US | #800-0198 | UK |
| #800-0197 | EU | #800-0199 | DEN |

- #508-0105 2" (5cm) Deluxe Pad
- #508-0152 2" (5cm) Softcare Pad



Pad sold separately



Radiolucent Armboard

- Slips under patient's shoulder and provides support during injection of contrast materials for MRI and CT scans or placement of PICC lines
- 2" conductive, vinyl-covered foam pad mounts with an attached sleeve and strap (sold separately)
- Radiolucent, Phenolic and Reusable

#800-0297 Armboard, Radiolucent, 6"W x 30"L

#508-0071 2" (5cm) Deluxe Pad

#508-0151 2" (5cm) Softcare Pad



Incline Armboard Pad

- Prevents hyperextension of the forearm
- Fits all 6"W x 26"L (15cm x 66cm) armboards
- Made with pressure-diffusing T-Foam™
- Covered in durable, black, conductive vinyl

#508-1355 Armboard, Incline Pad



Armboard Holder

- Stores two standard surgical table armboards
- Use convenient top shelf for clamps or other accessories
- Attaches to wall with four screws (supplied)

#800-0296 Armboard Holder



Schure Spine Frame

- Adjusts intraoperatively to the desired lordosis to get optimal surgical site exposure
- Transport dolly of durable stainless steel included
- Patient weight capacity is 500 lbs (227 kgs)
- Deluxe foam pads included

The Schure Spine Frame provides a stable and convenient platform for spine surgeries such as laminectomies, decompressions, and disc surgery. You can adjust it intraoperatively to the desired lordosis to get optimal surgical site exposure. Tapered adjustable pads provide optimal patient safety and comfort. The Schure Spine Frame has a patient weight capacity of 500 lbs (227 kgs) and is C-arm

compatible. It can be used with any OR table, and fits the Jackson Spine Table and Allen Flex Frame as well. Comes complete with transport storage dolly and crank handle, disposable head rest, and two arm cradles.

#800-0235 Schure Spine Frame

Optional Accessories

#508-0373 Disposable Laminectomy Arm Cradles, 12 prs/cs

#508-1340 Disposable Richard Slotted Head Rests, Extended, 6/cs

#508-0472 Schure Spine Frame Pad Set

#800-0277 Schure Spine Frame Crank Handle

#800-0265 Dolly, Schure Spine Frame



*TKR positions knee
in precise angles
and positions*



TKR Positioner for Total Knee Replacement

- *Track Locking Bracket securely locks the Base Plate to the table with two thumb screws*
- *Locks in any position with a simple turn of the wrist*
- *Smooth sliding motion for effortless leg positioning*
- *CamLoc Clamp sold separately (see below)*

Ball joint allows the surgeon to position the operative leg in extension, flexion, tilt and rotation.

- #800-0141** TKR Positioner, 24" (61cm)
- #800-0171** TKR Positioner, 30" (76cm)
- #800-0146** Disposable TKR Boot Liners, Sterile, 12/cs

Required Accessory: Choose One

CamLoc Clamp

- #800-0271** US **#800-0273** UK
- #800-0272** EU **#800-0274** DEN





Stainless Steel

Schure Foot for Total Knee Replacement

- This new design offers the same secure support as the stainless steel model without the concerns of radiopacity
- The surgeon can replace the knee while the patient's foot is held firmly in place
- Disposable conductive, non-latex foot pad included
- Now available in radiolucent, clear Lexan®
- Patient's well leg lies flat
- Schure Socket XL sold separately

- #800-0011 Schure Foot, Lexan
- #800-0062 Schure Foot, Stainless Steel
- #508-0111 Disposable Foot Pads, 12/cs
- #800-0134 Schure Socket XL



Lexan

Arthroscopic Leg Holder

- VELCRO® brand fastener strap wraps completely over leg holder for total stability
- Pressure relief foam pad included
- 5/8" (2cm) mounting post
- Schure Socket XL sold separately

Extremely versatile, easy-to-use and completely rigid. Spring-loaded braces slide into place and snap lock.

- #800-0023 Arthroscopic Leg Holder
- #508-0103 Disposable Premium Pads, 12/cs
- #800-0134 Schure Socket XL
- #508-0352 Well Leg Holder



Schure Grip Leg Holder

- Dual independent braces freely rotate 360° and spread horizontally from 3"-11"W (8cm-28cm) to fit any size leg at any angle
- Twist one brace side down, put the patient's leg on the leg holder, return the brace side back up, slide brace horizontally for a snug fit, and lock in place with convenient twist handle

Use also for ACL procedures. Stainless steel. 12"L x 1 1/2"D (31cm x 4cm) leg support x 8"H (20cm) braces, 16"H x 5/8"D (41cm x 1.6cm) mounting post. Schure Socket XL sold separately.

- #800-0065 Schure Grip Leg Holder
- #508-0103 Disposable Premium Pads, 12/cs
- #800-0134 Schure Socket XL





Telescoping mounting columns fit all tables

Foot Extension

- Mounts vertically as a foot support when using reverse Trendelenburg
 - Fits all OR tables and completely radiolucent
- Use as a table extension in the horizontal position to lengthen size of OR table. 2" (5cm) Deluxe Pad included.


#800-0036 Foot Extension, 20"W x 10"L (51cm x 25cm)

Required Simple Clamp: Order Two
#800-0228 US **#800-0264** UK
#800-0262 EU **#800-0263** DEN

Shoulder Supports

- Provides comfortable support for Trendelenburg positioning
- Deluxe pads on a broad support surface with ball and socket adjustments provide optimal positioning for any size patient
- Adjusts vertically to accommodate x-ray tops
- Schure Socket XL sold separately

#800-0100 Shoulder Supports
 Required Accessory
#800-0134 Schure Socket XL



Telescoping mounting columns fit all tables

Narrow Head Rest

- Improves surgeon's access at the head of the surgical table during ENT, neurosurgery, plastic surgery or ophthalmic surgery. Use instead of standard head rest.
- Radiolucent head rest includes 2" (5cm) premium pressure management pad. Deluxe pad included.

#800-0052 Narrow Head Rest
 Required Simple Clamp: Order Two
#800-0228 US **#800-0263** DEN
#800-0264 UK

Optional Accessory
#800-0086 Wrist Rest

Nissen Straps

- Offers support needed during steep reverse Trendelenburg positioning. Attaches to rail quickly and securely.
- 20"L x 5"W x 1" (51cm x 13cm x 3cm) thick pad offers complete support

#800-0055 Nissen Straps, 48"L x 2"W (19cm x 5 cm)
#800-0118 Bariatric Nissen Straps, 72"L x 2"W (183cm x 5cm)

Required Simple Clamp: Order Two
#800-0228 US **#800-0264** UK
#800-0262 EU **#800-0263** DEN




Table Width Extenders

- Attaches to OR side rail with mounting clamp
- Raise or lower height of extender to be level with the OR table, either with or without an x-ray top

Attached side rail allows any table accessory to be mounted. Deluxe pad included.

#800-0077 20"L x 4"W (51cm x 10cm)
#800-0079 15"L x 8"W (38cm x 20cm)
#800-0081 11"L x 8"W (28cm x 20cm)

Required Simple Clamp: Order One
#800-0228 US **#800-0264** UK
#800-0262 EU **#800-0263** DEN

Restraint Strap

- Amsco style
- All-rubber conductive strap
- 96"L (244cm) with two airplane-style buckles
- Very durable...easy to clean and disinfect
- Mounts to surgical table with or without hooks

#800-0057 Restraint Strap, No Hooks
#800-0058 Restraint Strap, 2 Hooks
#800-0119 Bariatric Restraint Strap, No Hooks
#800-0107 Bariatric Restraint Strap, w/Hooks





Multi-Axis Arm Positioner

- Versatile and comfortable positioning with the new ball socket locking system
- Arm platform measures 6"W x 16"L (15cm x 41cm)
- 2" (5cm) Deluxe Concave Foam Pad included
- 5/8" (2cm) Mounting Post

Perfect for seated, lateral, prone, neuro and park bench positions. Schure Socket XL sold separately.

- #800-0050** Multi-Axis Arm Positioner
- #800-0134** Schure Socket XL

Toboggan Arm/Leg Guard

- Available in four sizes: Small through X-Large
- Simply slide under OR table pad to secure arms or legs
- Physician can lean against guards without interfering with anesthesia lines

- #800-0084** Small, 6"H x 10"L (15cm x 25cm)
- #800-0083** Medium, 8"H x 16"L (20cm x 41cm)
- #800-0082** Large, 10½"H x 16"L (27cm x 41cm)
- #800-0276** X-Large, 10½"H x 24"L (27cm x 61cm)



Lateral Positioners

- Lateral Positioners are simple to set up and offer another alternative for lateral positioning
- Three swivel braces have Deluxe foam pads and rotate over the OR to conform to any size patient, they offer complete stability and ease of use
- 5/8" (2cm) stainless steel mounting post
- Schure Socket XL separately

- #800-0043** Lateral Positioners
- #800-0134** Schure Socket XL

Anesthesia Screen

- Adjusts from 15"-33" (38cm x 84cm) in height
- Horizontal adjustment feature facilitates viewing of patient
- Locks with a .25" (.6cm) twist of the handle
- Stainless steel construction
- Schure Socket XL sold separately

- #800-0016** Anesthesia Screen
- #800-0134** Schure Socket XL



Lateral Braces

- Secure and comfortable positioning for total hip replacement, and any other procedure requiring total lateral stability
- Three multi-axis braces easily accommodate any size patient for perfect lateral positioning
- Pressure relief foam with Deluxe cover provide extreme protection for the patient
- Schure Socket XL sold separately

- #800-0042** Lateral Braces
- #800-0134** Schure Socket XL

Leg Prepper

- Padded horseshoe cradle holds patient's leg comfortably during the prepping procedure
- Stainless steel 30"H x 5/8"D (76cm x 2cm) mounting post rotates, rises, and lowers easily to various prepping positions
- Schure Socket XL sold separately

- #800-0013** Leg Prepper
- #508-0107** Replacement Pad
- #800-0134** Schure Socket XL





Bi-Lateral Leg Prepper

- 500 lb (227 kg) patient weight capacity
- Pads remove for easy cleaning

The Bi-Lateral Leg Prepper safely and effectively holds both legs for bi-lateral preps. The unique rotating/pivoting design of the holders allow for infinite positioning off the extremity. Includes 2" (5cm) Deluxe Pads. Schure Socket XL sold separately.

#800-0236 Bi-Lateral Leg Prepper
#800-0134 Schure Socket XL



Well Leg Holder

- Helps prevent nerve damage
- Works on both left and right leg

The foam Well-Leg Holder easily supports with comfort the non-affected, well-leg during arthroscopic knee procedures, fracture repairs and cast procedures.

#508-0352 Well Leg Holder



Extremity Positioner

- Great for positioning the lower extremity for foot or ankle procedures
- Sealed seams for easy clean-up
- Dimensions: 16"L x 12"W x 7"H (41cm x 31cm x 18cm)

This pad is ideal and used regularly for prepping of the arm or leg surgical site.

#508-0353 Extremity Positioner

Universal Prepper

- Supports limbs of all sizes
- Pad removes for easy cleaning

The Universal Prepper safely and effectively holds upper or lower extremities. The unique pivoting design allows for infinite positioning of the extremity. Comes complete with 2" (5cm) thick Deluxe pad. Schure Socket XL sold separately.

#800-0161 Universal Prepper
#800-0134 Schure Socket XL



Saphenous Vein Harvest Pad

- Holds legs externally rotated and slightly flexed for minimally invasive, bridging or open harvesting procedures
- The cover is a multi-directional stretch cover. Inside consists of a top layer of conforming foam with a firm supportive base. Water-resistant cover is not made with natural rubber latex.
- Sealed watertight seams for easy clean-up

The Saphenous Vein Harvest Pad holds patient's legs at the optimal position for assessment and harvesting of the saphenous veins for coronary bypass surgery.

#508-0385 Saphenous Vein Harvest Pad

Elbow Arthroscopy Support

The Elbow Arthroscopy Support comfortably and securely holds the upper extremity for prepping and during surgical procedures. The unique pivoting design provides flexibility for all patient anatomies. Comes complete with Deluxe pad which can be removed for easy cleaning. Schure Socket XL sold separately.

#800-0168 Elbow Arthroscopy Support
#800-0134 Schure Socket XL





Irrigation Towers

- Can hold up to eight (8), five (5) liter bags
- Wide, inverted, 5-legged base resists tipping... no need for "sand-bagging"
- Adjustable from 5'- 8' (1.5m-2.4m) (9' with custom hangers)
- 2- or 4-channel designs
- Push-down squeeze release handles gives a 2-1 mechanical advantage
- No electric or hydraulic parts that require overnight recharging, heavy maintenance, or seals that leak

Great for fluid-intensive irrigation procedures such as arthroscopic, orthopedic, hysteroscopic, endoscopic, TURP, etc.

- #800-0157** Irrigation Tower, 4-Channel
- #800-0158** Irrigation Tower, 2-Channel



Safety Draw

- Avoid needle sticks
- Hands-free access to medication vials
- Easy viewing of medication labels
- Helps reduce medication errors

Safety-Draw maintains sterility while allowing access to non-sterile vials. This protects both patient and staff by helping to prevent dangerous medication mix-up, as well as needle sticks.

- #800-0059** Safety Draw, IV Pole Mount
- #800-0060** Safety Draw, Mayo Stand Mount

Schure Strap VELCRO® Brand Rolls

Schure Strap is the answer to your taping, strapping, tie-down and restraining requirements. It has VELCRO® brand fastener hook and loop back-to-back on 30' (9 meters) rolls and available in 4 widths. No sticky tape residue to damage equipment. Soft edge design to protect patient's skin. Not made with natural rubber latex.

- #800-0147** ¾"W (2cm)
- #800-0148** 1½"W (4cm)
- #800-0149** 3"W (8cm)
- #800-0150** 5"W (13cm)



A-Line Positioner w/Straps

A-Line Positioner w/Straps consist of two exclusive 1"W x 17"L (3cm x 43cm) soft edge straps which reduce the risk against skin ulceration. This anatomically formed base provides the best angle while keeping the skin taut for optimal needle insertion.

- #800-0143** Disposable A-Line Positioner w/Straps, 12/cs



A-Line Positioners

- Stable and accurate positioning
- Enhanced comfort for the patient
- No need to try to find a make-shift positioner
- Has slot to secure to armboard

Disposable support positioner to provide for a fast, consistent and reliable A-Line insertion.

- #800-0144** Disposable A-Line Positioners, 30/cs





Schure Socket XL

- Hardened clamping surface creates maximum clamping force stopping movement of all accessories
- Attaches along the side rail, without the need of notches

Accepts any round mounting posts up to 5/8"D (1.6cm) and rectangular mounting posts up to 1" x 3/8" (2.5cm x .95cm). Quickly and easily connects and disconnects securely with a quick turn anywhere along the side rail, even over side rail stand off.

Schure Socket XL

#800-0134 US #800-0200 UK
#800-0244 EU #800-0201 DEN

800 lb. (363 kg)
Patient Wt. Capacity

Patented

Schure Socket

- Fits round mounting posts up to 5/8"D (2cm)
- Locks securely with a twist of the ergonomic handle
- Attaches along the side rail, without the need of notches

The industrial strength Schure Socket fits all OR tables with standard 3/8"W x 1 1/8"H (.95cm x 2.85cm) side rails.

#800-0006 US #800-0203 UK
#800-0202 EU #800-0204 DEN



800 lb. (363 kg)
Patient Wt. Capacity

Patented



SpringLoc Clamp

- Securely locks accessories with a flat 1" x 1/4" (2.5cm x .6cm) mounting post
- Attaches along the side rail, without the need of notches
- Quick & easy pull-to-release lever

#800-0127 SpringLoc Clamp

800 lb. (363 kg)
Patient Wt. Capacity

Universal Rail Clamp

- Spring-loaded and attaches to OR table side rail anywhere along the rail
- Accessories can be mounted vertically or horizontally
- Locks securely with a twist of the ergonomic handle

Accepts both 5/8"D (1.6cm) and 1" x 3/8" (2.5cm x .95cm) flat bar mounting posts.

#800-0085 Universal Rail Clamp



350 lb. (158 kg)
Patient Wt. Capacity

Patented

Simple Clamp

- Attaches at OR table rail ends or notches
- Accessories can be mounted vertically or horizontally

Securely locks all accessories with flat 1" x 3/8" (2.5cm x .95cm) vertically mounting blades anywhere along the table side rail. 1" x 1/4" (2.5cm x .63cm) horizontally mounting blades.

#800-0228 US #800-0264 UK
#800-0262 EU #800-0263 DEN



720 lb. (327 kg)
Patient Wt. Capacity



CamLoc Clamp

- It can attach to the rail at any location over multiple drapes and accepts accessories with round posts up to 5/8"D (1.6cm) and flat posts 1" x 3/8" (2.5cm x .95cm)
- Attaches along the side rail, without the need of notches

The CamLoc Clamp is designed for use in the sterile field.

CamLoc Clamp

#800-0271 US #800-0273 UK
#800-0272 EU #800-0274 DEN


500 lb. (227 kg)
Patient Wt. Capacity

Deluxe Rail Clamp

- Attaches at OR table rail ends or notches
- Accessories can be mounted vertically or horizontally

Securely locks all accessories with 1" x 3/8" (2.5cm x .95cm) flat mounting post and round mounting posts up to 5/8"D (1.6cm) mounting blades anywhere along the table side rail.

#800-0248 US #800-0206 UK
#800-0205 EU #800-0207 DEN



500 lb. (227 kg)
Patient Wt. Capacity



Comfort Mat


The 20" x 33" (51cm x 84cm) Comfort Mat is ideal at scrub sinks, nurse's stations or at OR tables. Helps reduce fatigue, pain and discomfort associated with standing for long periods of time.

- Absorbent-resistant closed-cell foam
- Beveled, low-profile edges
- Made in USA
- Non-slip surface is chemical resistant
- Custom sizes are available

#800-0145 Comfort Mat



Schure Slide XLT



Schure Slide XL



Schure Slide, Ergo



Schure Slide, Anti-Static



Schure Slide, Full-Size

SCHURE SLIDE
patient transfer system

Schure Slide is smooth surfaced, flexible, and extremely durable. After being placed under a patient, personnel can effortlessly slide the patient from one gurney or table to another. The Full-Size Schure Slide enables personnel to easily slide the patient from one bed to another. Choose Flexible or Rigid. Patented.


- Reduce back injuries and increase comfort
- Non-disposable, no need to waste red bags
- Won't stain—spray and wipe clean
- Speed up turnover time—hang Schure Slide in every OR
- Radiolucent & not made w/natural rubber latex

ORDERING INFORMATION

	1+	5+	10+
#800-0072 Schure Slide, XL, 33"L x 22"W x 1/16" (84cm x 56cm x .2cm)			
#800-0071 Schure Slide, Rigid, 33"L x 22"W x 3/16" (84cm x 56cm x .5cm)			
#800-0162 Schure Slide, Ergo, 33"L x 23"W x 1/16" (84cm x 58cm x .2cm)			
#800-0067 Schure Slide, Anti-Static, Black, 33"L x 22"W x 1/16" (84cm x 56cm x .2cm)			
#800-0068 Schure Slide, Anti-Static Rigid, Black, 33"L x 22"W x 3/16" (84cm x 60cm x .5cm)			
#800-0069 Schure Slide, Full Size, 72"L x 22"W x 3/16" (183cm x 56cm x .5cm)			
#800-0070 Schure Slide, Full Size, Anti-Static, Black, 72"L x 22"W x 3/16" (183cm x 56cm x .5cm)			
#800-0120 Schure Slide, XLT, 72"L x 22"W x 3/16" (183cm x 56cm x .5cm)			



Reduces
Chance of Cross-
Contamination
Slip-Resistant




Ordering Information

Reusable antimicrobial coated positioners ideal for extremity positioning. The positioners are radiolucent and the coating is impregnated with Agion antimicrobial ionic silver that inhibits the growth of bacteria.

Item #	Description	Dimensions
A. #508-0355	Double Angle Wedge, 35°/55°	24"L x 9¾"W x 6¾"H (61cm x 25cm x 17cm)
#508-0356	Angle Wedge, 45°	24"L x 7"W x 7"H (61cm x 18cm x 18cm)
#508-0357	Wedge, 15°	10½"L x 9¾"W x 2¾"H (27cm x 25cm x 7cm)
B. #508-0358	Double Angle, 55°/35°	14"L x 9½"W x 6½"H (36cm x 24cm x 17cm)
#508-0359	Wedge, 20°	11"L x 8½"W x 3¾"H (28cm x 22cm x 10cm)
#508-0360	Head & Neck	10¾"L x 7"W x 6"H (27cm x 18cm x 15cm)
C. #508-0361	Head Support, Occipital Adult	10"L x 10"W x 4"H (25cm x 25cm x 10cm)
D. #508-0364	Semicircle	14"L x 8½"W x 4½"H (36cm x 22cm x 11cm)
#508-0365	Rolls, Small	24"L x 6"D (61cm x 15cm)
#508-0366	Rolls, Large	24"L x 8"D (61cm x 20cm)
E. #508-0367	Donut, Large, 6 pk.	9"D x 1½"H (23cm x 4cm)
#508-0368	Donut, Small, 6 pk.	6"D x 1½"H (15cm x 4cm)
F. #508-0362	Rectangular, Large	14"L x 9¼"W x 2"H (36cm x 24cm x 5cm)
#508-0363	Rectangular, Small	9½"L x 7"W x 2"H (24cm x 18cm x 5cm)



Variety of Shapes
& Sizes Available
High-Density
Polyfoam



Ordering Information

These high-density Polyfoam positioning pads are supportive, comfortable and economical for all your positioning requirements. Available in a variety of shapes and sizes. Custom sizes available. Please call for a quote.

Item #	Description	Dimensions
A #508-0369	Disposable, Wedge	24"L x 9¾"W x 6¾"H (61cm x 25cm x 17cm)
#508-0386	Disposable, Wedge	24"L x 7"W x 7"H (61cm x 18cm x 18cm)
#508-0387	Disposable, Body Aligner, Wedge	28"L x 8"W x 10"H (71cm x 20cm x 25cm)
#508-0388	Disposable, Body Aligner, Half Wedge	14"L x 8"W x 10"H (36cm x 20cm x 25cm)
B #508-0379	Disposable, Chest Roll, 12/cs	5"D x 17"L (13cm x 43cm)
C #508-0381	Disposable, Adult Head Support, 12/cs	9"L x 8"W x 3½"H (23cm x 20cm x 9cm)
D #508-0374	Disposable, Richards Head Rest, 6/cs	8½"L x 8"W x 3¾"H (22cm x 20cm x 10cm)
E #508-0375	Disposable, Richards Slotted Head Rest, 6/cs	8½"L x 8"W x 3¾"H (22cm x 20cm x 10cm)
F #508-0382	Disposable, Pediatric Head Support, 24/cs	6½"L x 7"W x 3¼"H (17cm x 18cm x 8cm)
G #508-0380	Disposable, Convuluted Sheet, 12/cs	20"L x 12"W x 2"H (51cm x 31cm x 5cm)
H #508-0373	Disposable, Laminectomy Arm Cradles, 12 pr/cs	5"W x 24"L x 3"H (13cm x 61cm x 8cm)
I #508-0371	Disposable, Small Donut, 36/cs	1½"H x 6"D (4cm x 15cm)
J #508-0372	Disposable, Super Donut, 36/cs	9"W x 10"L x 1"H (23cm x 25cm x 3cm)
K #508-0370	Disposable, Large Donut, 24/cs	1½"H x 9"D (4cm x 23cm)



Reusable, E-Z clean patient positioners

Gel Head Positioners

- 1 Gel Donut**
#508-0431 Adult
 3" I.D x 8" O.D x 1 3/4"H (8cm x 20cm x 4cm)
#508-0432 Tall Adult
 4" I.D x 7 3/4" O.D x 3"H (8cm x 20cm x 8cm)
#508-0446 Pediatric
 2 1/4" I.D x 5 1/2" O.D x 1 1/4"H (8cm x 20cm x 8cm)
#508-0478 Neonatal
 1 1/2" I.D x 3 1/4" O.D x 3/4"H (4cm x 9cm x 2cm)
- 2 Gel Horseshoe**
#508-0443 Adult
 3" I.D x 8" O.D x 1 3/4"H (8cm x 20cm x 3cm)
#508-0444 Tall Adult
 4" I.D x 7 3/4" O.D x 3"H (8cm x 20cm x 8cm)
#508-0447 Pediatric, Medium
 2 1/4" I.D x 5 1/2" O.D x 1 1/4"H (6cm x 14cm x 3cm)
#508-0479 Neonatal
 1 1/4" I.D x 3 1/4" O.D x 3/4"H (4cm x 9cm x 2cm)
- 3 Gel Head Pillow w/Centering Dish**
#508-0470 Lateral
 10"W x 9"L x 1"H (25cm x 23cm x 3cm)
#508-0473 Supine
 10"W x 9"L x 1 3/4"H (25cm x 23cm x 4cm)
#508-0486 Pediatric Supine
 7" x 6" x 2" (18cm x 15cm x 5cm)

- #508-0487** Pediatric Lateral
 7" x 6" x 1" (18cm x 15cm x 3cm)
#508-0488 Neonatal Supine
 5" x 4" x 2" (13cm x 10cm x 5cm)
#508-0489 Neonatal Lateral
 5" x 4" x 1" (13cm x 10cm x 3cm)
- Gel Head Rest**
#508-0520 Contoured
 7" x 7" x 3/4" (18cm x 18cm x 8cm)
#508-0532 Ophthalmic
 12" x 10" x 3/8" (31cm x 25cm x 10cm)
- 4 Gel Richard Type Headrest**
#508-0462 T Shape
 9" x 8" x 4" (23cm x 20cm x 10cm)
#508-0463 Tall T Shape
 9" x 8" x 4 1/2" (23cm x 20cm x 11cm)
- 5 Gel Prone/Supine Headrest**
#508-0464 w/Anesthesia Tube Cutouts
 11" x 9" x 6" (28cm x 23cm x 15cm)

Chest and Axillary Positioners & Rolls

- 6 Gel Pediatric Chest Roll Flat Bottom**
#508-0418 10" x 2" x 2" (25cm x 5cm x 5cm)
#508-0419 10" x 3" x 3" (25cm x 8cm x 8cm)
#508-0420 4" x 2" x 2" (10cm x 5cm x 5cm)
#508-0421 4" x 3" x 3" (10cm x 8cm x 8cm)
#508-0422 Small
 6 1/2"L x 2 1/2"W x 2"H (17cm x 6cm x 5cm)
#508-0423 6" x 3" x 3" (15cm x 8cm x 8cm)




Gel Pads & Positioners to position every patient

- #508-0424** Small
 8"L x 2"W x 2"H (20cm x 5cm x 5cm)
#508-0425 8" x 3" x 3" (20cm x 8cm x 8cm)
#508-0550 12" x 2" x 2" (31cm x 5cm x 5cm)
#508-0551 12" x 3" x 3" (31cm x 8cm x 8cm)
- Gel Chest Roll-Flat Bottom**
#508-0456 12"L x 4"W x 3 3/4"H (31cm x 10cm x 10cm)
#508-0457 12" x 6" x 6" (31cm x 15cm x 15cm)
#508-0458 16"L x 4"W x 3 3/4"H (41cm x 10cm x 10cm)
#508-0459 16" x 6" x 6" (41cm x 15cm x 15cm)
#508-0460 20"L x 4"W x 3 3/4"H (51cm x 10cm x 10cm)
#508-0461 20"L x 7"W x 5 1/2"H (51cm x 18cm x 15cm)
- Gel Contoured Chest Roll Flat Bottom**
#508-0529 20" x 6 1/4" x 5" (51cm x 16cm x 13cm)
#508-0530 16" x 6 1/4" x 5" (41cm x 16cm x 13cm)
#508-0534 12" x 6 1/4" x 5" (31cm x 16cm x 13cm)
- 8 Gel Chest Positioner**
#508-0474 Left, 2 3/28"L x 8"W x 2"H (2cm x 20cm x 5cm)
#508-0475 Right, 2 3/28"L x 8"W x 2"H (2cm x 20cm x 5cm)

- #508-0476** Left & Right, 2 3/28"L x 8"W x 2"H (2cm x 20cm x 5cm)
Gel Lateral Positioner Pad
#508-0494 28" x 20" x 1/2" (71cm x 51cm x 1cm)
- 10 Gel Axillary Roll Cover w/Hook & Loop VELCRO® brand fastener**
#508-0556 17" x 2" x 1/2" (43cm x 5cm x 1cm)
#508-0557 18" x 18" x 1/2" (46cm x 46cm x 1cm)


Foot, Hip, Heel, Knee & Stirrup Pads

- 11 Gel Candy Cane Stirrup Sheaths**
#508-0453 24" L x 2" O.D (61cm x 5cm)
- 12 Gel Ankle Candy Cane Stirrup Straps**
#508-0506 Small, 12"L x 1"W x 1/2" (32cm x 5cm x 1cm)
#508-0507 Large, 17 3/8" x 1 7/8" x 3/8" (32cm x 5cm x 1cm)
- 13 Gel Clam Shell Heel Pad w/Velcro®**
#508-0451 12 1/2" x 7 3/4" x 1/2" Thick (32cm x 20cm x 1cm)
- 14 Gel Lithotomy Boot Pad**
#508-0437 Small, 24" x 11 1/2" x 6" (61cm x 29cm x 15cm)
#508-0438 Large, 27" x 14 1/2" x 8" (69cm x 37cm x 20cm)



Maximum protection against surgical patient pressure ulcers

Our patient positioners lower surgical risk safely and securely




Foot, Hip, Heel, Knee & Stirrup Pads

- 15 Gel Knee Crutch Pad**
- #508-0452** Small, 11"W x 9"L x 3/8"H (28cm x 23cm x 1cm)
- #508-0415** Med., 17"W x 11"L x 1/4"H (43cm x 28cm x 1cm)
- #508-0435** Lg., 20"L x 20"W x 1/2" (51cm x 51cm x 1cm)
- #508-0465** X-Lg., 27" x 11" x 1/2" (69cm x 28cm x 1cm)
- Gel Bariatric Knee Crutch Pad**
- #508-0483** Small, 11" x 9" x 1" (28cm x 23cm x 3cm)
- #508-0484** Medium, 17" x 11" x 1" (43cm x 28cm x 3cm)
- #508-0485** Large, 20" x 20" x 1" (51cm x 51cm x 3cm)
- 16 Gel Foot/Heel/Sole Protector Pad w/Velcro®**
- #508-0450** 6"W x 15"L x 1/4"H (38cm x 10cm x 1cm)
- 17 Orthopedic Hip Positioner with Rigid Foam Core**
- #508-0561** 16" x 8 1/2" x 2 5/8" (41cm x 22cm x 7cm)
- 18 Gel Sandbags**
- #508-0454** Standard 10"L x 4 1/2"W x 2"H (25cm x 11cm x 5cm)

- #508-0455** Large 14"L x 4 1/2"W x 2 1/4"H (36cm x 11cm x 6cm)
- #508-0466** X-Large 16" x 12" x 1 1/2" (41cm x 31cm x 4cm)
- #508-0496** Pediatric 6" x 1 1/2" x 1 1/2" (15cm x 4cm x 4cm)
- #508-0497** Pediatric 8" x 1 1/2" x 1 1/2" (20cm x 4cm x 4cm)
- #508-0498** Pediatric 6" x 2" x 2" (15cm x 5cm x 5cm)
- #508-0499** Medium 12" x 4 1/2" x 1 1/2" (31cm x 11cm x 4cm)
- #508-0500** Medium Extra Thick 12" x 4 1/2" x 2" (31cm x 11cm x 5cm)
- #508-0511** Large 18" x 4 1/2" x 2" (46cm x 11cm x 5cm)
- #508-0519** X-Large 21" x 5" x 2 1/2" (53cm x 13cm x 6cm)
- #508-0521** Neonatal 4" x 1 1/2" x 1 1/2" (10cm x 1cm x 1cm)
- #508-0522** Neonatal 6" x 1 1/2" x 1 1/2" (15cm x 1cm x 1cm)
- #508-0523** Neonatal 8" x 1 1/2" x 1 1/2" (20cm x 1cm x 1cm)
- #508-0524** Neonatal 4" x 1" x 1" (10cm x 3cm x 3cm)
- #508-0525** Neonatal 6" x 1" x 1" (15cm x 3cm x 3cm)
- #508-0526** Neonatal 8" x 1" x 1" (20cm x 3cm x 3cm)

- #508-0528** Pediatric 12" x 1 1/2" x 1 1/2" (31cm x 4cm x 4cm)
- #508-0547** Medium, Bariatric 13" x 7" x 3" (33cm x 18cm x 8cm)
- #508-0548** Large, Bariatric 15" x 8" x 3" (38cm x 20cm x 8cm)
- Elbow, Arm & Shoulder Pads**
- 19 Heel Cup/Mini Positioner**
- #508-0477** 6"L x 3"W x 1 1/2"H (15cm x 8cm x 4cm)
- 20 Oval Ulnar Nerve Protector w/Velcro®**
- #508-0449** 6"W x 13"L x 1 1/2"H (15cm x 33cm x 1cm)
- 21 Gel Square Ulnar Nerve Pad**
- #508-0513** 24" x 21" x 1/4" (61cm x 53cm .6cm)
- #508-0518** Gel Ulnar Nerve Protector with VELCRO® brand fastener, 30" x 25" x 1/4" (76cm x 64cm x .6cm)
- #508-0448** Gel Square Ulnar Nerve Pad 17" x 11" x 3/8" (43cm x 28cm x 1cm)
- 22 Gel Armboard Pad**
- #508-0433** Large, Full 24"L x 4 1/2"W x 1/4"H (61cm x 11cm x 1cm)
- #508-0434** X-Large 26" x 6" x 3/8" (66cm x 15cm x 1cm)
- #508-0445** Small 15" x 4" x 1/4" (38cm x 10cm x .6cm)

- #508-0514** Large 25" x 7" x 1" (64cm x 18cm x 3cm)
- #508-0535** Medium 20" x 5" x 1/4" (51cm x 13cm x .6cm)
- #508-0480** Gel DuraBoard Pad Short 15" x 6" x 3/8" (38cm x 15cm x 1cm)
- 23 Gel Contoured Armboard Pad**
- #508-0508** Medium 20"L x 5 3/4"W x 2 1/4" H (51cm x 15cm x 6cm)
- #508-0512** Small 15"L x 5 1/2"W x 2"H (38cm x 14cm x 5cm)
- 24 Gel Shoulder/Kidney Brace Pad Set**
- #508-0509** 7" x 6" x 1" (18cm x 15cm x 3cm)
- 25 Gel Elbow Pad**
- #508-0555** 3 1/4" O.D. X 1 1/2" X 3/4" (8cm x 4cm x 2cm)
- Surgical Table Pads**
- 26 Gel Three Piece Pad Set Universal**
- 27 Gel Head Section Pad**
- #508-0427** 20" x 9" x 1/2" (51cm x 23cm x 1cm)
- 28 Gel Torso Section Pad**
- #508-0428** 39 1/2" x 20" x 1/2" (100cm x 51cm x 1cm)




(Surgical Table Pads continued)

- 29**
#508-0429 Gel Foot Section Pad w/Cutout
20" x 20" x 1/2" (51cm x 51cm x 1cm)
- 30**
#508-0442 Gel Full Table Pad
72"L x 20"W x 1/2"H (183cm x 51cm x 1cm)
- 31**
#508-0440 Gel 3/4 Table Pad
46"L x 20"W x 1/2"H (117cm x 51cm x 1cm)
- #508-0536** Gel Cysto Table Pad
44" x 28" x 1/2" (112cm x 71cm x 1cm)
- #508-0545** Gel Cysto Table Pad
30" x 24" x 1/2" (76cm x 61cm x 1cm)
- 32**
#508-0441 3/4 Table Pad w/GU Cutout
46"L x 20"W x 1/2"H (117cm x 51cm x 1cm)
- 33**
#508-0517 Gel Cysto Table Pad
25 1/2" x 11 1/2" x 1/2" (65cm x 29cm x 1cm)
- #508-0533** Gel Table Pad
18 1/2" x 7" x 2 1/2" (47cm x 18cm x 6cm)
- #508-0542** Gel Table Pad
36" x 10" x 1/2" (91cm x 25cm x 1cm)
- #508-0543** Gel Table Pad
60" x 10" x 1/2" (152cm x 25cm x 1cm)
- #508-0544** Gel Flat Pad
33" x 18" x 1/2" (84cm x 46cm x 1cm)
- #508-0546** Gel Table Pad
48" x 20" x 1" (122cm x 51cm x 3cm)

- 34**
#508-0537 Gel Fracture Hip Pad
9" (23cm) tapers to 3" x 7 x 1/2" w/1 3/4"
(8cm x 18cm x 1cm w/4cm) Post Hole
- 35**
#508-0539 Gel Sacral Pad
9"(23cm) tapers to 3" x 7 x 1/2"
(8cm x 18cm x 1cm)
- 36**
#508-0538 Gel Fracture Table Post Pad w/Hook &
Loop VELCRO® brand fastener,
7" x 5 1/2" x 1/4" (18cm x 14cm x .6cm)

Specialty Pads

- 37** **Gel Wilson Frame Pad**
- #508-0471** Each
29" x 10" x 1/2" (74cm x 25cm x 1cm)
- #508-0558** With Velcro®, Set of 2
29"L x 10"W x 1/2"H (74cm x 25cm x 1cm)
- 38**
#508-0430 Gel Hip Pad
20" x 20" x 1/2" (51cm x 51cm x 1cm)
- #508-0502** Gel Table Peg Board Pad
32" x 20" x 1/2" (81cm x 51cm x 1cm)
- 39**
#508-0436 Gel Pediatric Pad
17" x 11" x 1/4" (43cm x 28cm x .6cm)
- #508-0503** Gel Pad
25" x 12" x 3/4" (64cm x 31cm 2cm)
- #508-0469** Gel Layer Pad
16" x 8" x 1/2" (41cm x 20cm x 1cm)
- #508-0527** Gel Pad
30" x 10" x 1/2" (76cm x 25cm x 1cm)



Maximum protection
against surgical
patient pressure ulcers

Our patient positioners
lower surgical risk
safely and securely

- 40**
#508-0504 Trapezoidal Knee Positioner
10"L x 7"W x 3 1/2"H (25cm x 18cm x 9cm)
- 41**
#508-0491 Gel Graduated Thoracic Positioner
14" x 9" x 2 1/2"-1 1/2"
(36cm x 23cm x 6cm-4cm)
- 42**
#508-0510 Gel Arm Toboggan Pad
7 3/4" x 6 1/4" x 3/8" (20cm x 16cm x 1cm)
- #508-0467** Gel Leg Toboggan Pad
17" x 6 1/2" x 3/8" (43cm x 17cm x 1cm)
- #508-0553** Gel Pad for Arm Shield w/pocket att.
20" L (51cm)
- #508-0554** Gel Toboggan Pad Shield
16 1/2" x 7 1/2" x 3/8" (42cm x 19cm x 1cm)
- 43**
#508-0490 Gel Andrews Spinal Frame Pad
22" x 16" x 1" (56cm x 41cm 3cm)
- 44**
#508-0492 Gel Neonatal Bassinet Pad
26" x 12" x 1/2" (66cm x 31cm x 1cm)
- #508-0493** Gel Neonatal Isolette Pad
26" x 6" x 1/2" (66cm x 15cm x 1cm)

- 45**
#508-0501 Gel Table Bean Bag Pad
36" x 28" x 1/2" (91cm x 71cm x 1cm)
- 46**
#508-0505 Gel Mayo Stand Cover
21" x 18 1/2" x 5/8" (53cm x 47cm x 2cm)

Miscellaneous Pads (no Pictures)

- #508-0540** Gel Commode Seat Pad w/Cutout
18" x 18" x 1" (46cm x 46cm x 3cm)
- #508-0541** Gel Commode Seat Pad
16" x 16" x 5/8" (41cm x 41cm x 2cm)
- #508-0549** Gel Summersville Pad
12" x 12" x 1/4" (31cm x 31cm x .6cm)
- #512-0047** Gel Repair Patch
3" x 6" (8cm x 15cm)
- #508-0439** Prone and Lateral Positioner
15" x 7" x 1/4" (38cm x 18cm x .6cm)
- #508-0481** Med. Bariatric Lap Belt Pad
36" x 6" x 1/4" (91cm x 15cm x .6cm)
- #508-0482** Large Bariatric Lap Belt Pad
48" x 7" x 1/4" (121cm x 18cm x .6cm)
- #508-0495** Epidural Head Rest
13" x 5" x 1/4" (33cm x 13cm x .6cm)
- #508-0531** OR Table Pad
37" x 12" x 1" (127cm x 31cm x 2cm)


Surfaces & Pads

SchureMed has over 20 years experience in the area of pressure management. We have worked with engineered foams and covering materials to provide you with the best combinations for your pressure management requirements in the surgical suite and surrounding areas. We realize that pressure sores are getting a lot of attention by healthcare professionals and SchureMed is your best option. No single foam offers effective pressure management, therefore SchureMed offers a variety of options to meet your specific needs and budget.

Research findings... why is it important to my OR?

PREVALENCE


- One research study concluded that 23% of all nosocomial pressure ulcers are related to the OR experience¹
- In an experiment involving 323 patients between 2000 and 2004, 38% of patients who had the standard pad during surgery developed the sores. But only 7% of patients who had the pressure-relieving pad developed pressure ulcers.³



Hospital Acquired Pressure Ulcers

COST TO TREAT

- The annual cost to treat surgical patients' pressure ulcers — nearly 42% of all hospital-acquired pressure ulcers occur in surgical patients — is between \$750 million and \$1.5 billion.²
- The treatment of a single ulcer has been estimated to cost \$14,000 to \$40,000.⁶



Annual Cost of Hospital—Acquired Ulcers

WHO IS AT RISK? AORN RECOMMENDATION

- The risk of intraoperative ulcerations increases as surgical time increases... Therefore all surgical patients undergoing prolonged procedures should be considered at risk for intraoperative ulceration.⁴

All surgical patients should be considered at-risk for pressure ulcer development; therefore, perioperative departments should develop and implement strategic plans for pressure ulcer prevention.⁵

SCHUREMED PAD TECHNOLOGY

- Available with glued or welded seams for easy cleaning and disinfection
- Available in OR Table Pads, Stretcher Pads and Positioners
- Comes complete with VELCRO® brand fastener Kit when applicable
- Available in 2"-5" (5cm-13cm) thicknesses
- Not made with natural rubber latex

¹ Beckrich K., Aronovitch SA Hospital- acquired pressure ulcers: a comparison of costs in medical vs. surgical patients. *Nurs Econ.* 1999; 17(5):263-71

² Beckrich K, Aronovitch SA. Hospital-acquired pressure ulcers: a comparison of costs in medical versus surgical patients. *Nurs Econ.* 1999;17(5):263-271.


³ Scott-Williams S. Perioperative Pressure Ulcer Prevention Program (PPUPP): An innovative effort to prevent pressure ulcers in surgical patients: American Academy of Nursing. http://www.aannet.org/files/public/PPUPP_template.pdf

⁴ Aronovich SA Intraoperatively acquired pressure ulcer prevalence: a national study. *J Wound Ostomy. Conference Nurs.* 1999;26(3) 130-136


⁵ AORN J 89 (March 2009) 538-548. © AORN, Inc, 2009.

⁶ The American Hospital Association. Hospital Statistics. Chicago, Ill: Healthcare Infosource, Inc;1998.21. Barczak C, Barnett RI, Childs EJ, Bosley LM. Fourth pressure ulcer prevalence survey. *Adv Wound Care.* 1997;10(4):18-26.

2" Standard Pad vs 2" Standard with Gel Overlay vs 2" Softcare vs 2" Softcare with Gel



Pressure Maps: Perioperative Surfaces



Subject: 5' 9" Male weighing 200 lbs (91 kgs) in the supine position



Standard

Consists of standard open cell high density polyurethane foam with a 3-Ply PVC Conductive Cover.



Deluxe

Consists of standard open cell high density polyurethane foam with a multidirectional stretch anti-shear 4-Way Stretch Conductive Cover.



Deluxe Gel

Consists of standard open cell high density polyurethane foam support layer with a top layer of polymer gel all covered with a multidirectional stretch anti-shear 4-Way Stretch Conductive Cover.



Softcare

Consists of high density polyurethane foam support layer with a top layer of memory foam covered with a multidirectional stretch anti-shear 4-Way Stretch Conductive Cover.



Softcare Plus

Consists of a closed cell rubber anti-bottoming layer then a high density polyurethane foam support layer with a top layer of memory foam covered with a multidirectional stretch anti-shear 4-Way Stretch Conductive Cover.



Softcare Silver

Consists of a closed cell rubber anti-bottoming layer then a high density polyurethane foam support layer with a top layer of memory foam covered with a multidirectional stretch anti-shear 4-Way Stretch Conductive Cover with embedded silver which inhibits growth of odor-causing bacteria.




Softcare Plus Silver


Consists of a closed cell rubber anti-bottoming layer then a high density polyurethane foam support layer with a top layer of memory foam covered with a multidirectional stretch anti-shear 4-Way Stretch Conductive Cover with embedded silver which inhibits growth of odor-causing bacteria.




A. Berchtold B810, B710




I. Midmark Chic CST




Q. Skytron 3501B




B. Berchtold Operon D850, D750




J. Midmark Montreal




R. Skytron 3600




C. Maquet/Siemens 1130




K. Getinge/Shampaine 1400/1500/1900




S. Skytron 6700, 6701




D. Maquet 1130 Split Leg




L. Getinge/Shampaine 1501 Radi-Op




T. Steris-Amsco 1080, 2080




E. Maquet 1130 Orthopedic




M. Getinge/Shampaine 4900/5100B, 5100E




U. Steris-Amsco 1080, 2080 w/Leg Ext.




F. Maquet Alphamaxx



N. Skytron 5001, 6001, 6500, 6600




V. Steris-Amsco 3080, 3085




G. Maquet Alphastar 1132




O. Skytron 3100




W. Steris-C-Max




H. Midmark 7100



P. Skytron 3500



X. Steris Surgi-Graphic 6000



		STANDARD		DELUXE		DELUXE GEL				
		US	2"	3"	2"	3"	2.25"	3.25"	4.25"	5.25"
		METRIC	5cm	8cm	5cm	8cm	6cm	8cm	11cm	13cm
	PRICE									
A	Berchtold B810, B710	#508-0169	#508-0170	#508-0171	#508-0172	#508-0189	#508-0190	#508-0191	#508-0192	
B	Berchtold Operon D850, D750	#508-0193	#508-0194	#508-0195	#508-0196	#508-0213	#508-0214	#508-0215	#508-0216	
C	Maquet/Siemens 1130	#508-0016	#508-0217	#508-0017	#508-0218	#508-0342	#508-0233	#508-0234	#508-0235	
D	Maquet 1130 Split Leg	#508-0236	#508-0237	#508-0238	#508-0239	#508-0256	#508-0257	#508-0258	#508-0259	
E	Maquet 1130 Orthopedic	#508-0260	#508-0261	#508-0262	#508-0263	#508-0280	#508-0281	#508-0282	#508-0283	
F	Maquet Alphamaxx	#508-0284	#508-0285	#508-0286	#508-0287	#508-0304	#508-0305	#508-0306	#508-0307	
G	Maquet Alphastar 1132	#508-0308	#508-0309	#508-0310	#508-0311	#508-0328	#508-0563	#508-0564	#508-0565	
H	Midmark 7100	#508-0013	#508-0566	#508-0014	#508-0567	#508-0341	#508-0582	#508-0583	#508-0584	
I	Midmark Chic CST	#508-0585	#508-0586	#508-0587	#508-0588	#508-0605	#508-0606	#508-0607	#508-0608	
J	Midmark Montreal	#508-0609	#508-0610	#508-0611	#508-0612	#508-0629	#508-0630	#508-0631	#508-0632	
K	Getinge/Shampaine 1400/1500/1900	#508-0034	#508-1297	#508-0035	#508-1298	#506-0348	#508-1313	#508-1314	#508-1315	
L	Getinge/Shampaine 1501 Radi-Op	#508-0028	#508-0633	#508-0029	#508-0634	#508-0346	#508-0649	#508-0650	#508-0651	
M	Getinge/Shampaine 4900/5100B, 5100E	#508-0031	#508-0652	#508-0032	#508-0653	#508-0347	#508-0668	#508-0669	#508-0670	
N	Skytron 5001, 6001, 6500, 6600	#508-0025	#508-0690	#508-0026	#508-0691	#508-0345	#508-0706	#508-0707	#508-0708	
O	Skytron 3100	#508-0019	#508-0709	#508-0020	#508-0710	#508-0343	#508-0725	#508-0726	#508-0727	
P	Skytron 3500	#508-0022	#508-0728	#508-0023	#508-0729	#508-0344	#508-0744	#508-0745	#508-0746	
Q	Skytron 3501B	#508-0395	#508-0747	#508-0396	#508-0748	#508-0399	#508-0763	#508-0764	#508-0765	
R	Skytron 3600	#508-0766	#508-0767	#508-0768	#508-0769	#508-0786	#508-0787	#508-0788	#508-0789	
S	Skytron 6700, 6701	#508-0390	#508-0790	#508-0391	#508-0791	#508-0394	#508-0806	#508-0807	#508-0808	
T	Steris-Amsco 1080, 2080	#508-0007	#508-0809	#508-0008	#508-0810	#508-0339	#508-0825	#508-0826	#508-0827	
U	Steris-AMSCO 1080/2080 w/Leg Ext.	#508-0010	#508-0828	#508-0011	#508-0829	#508-0340	#508-0844	#508-0845	#508-0846	
V	Steris-Amsco 3080, 3085	#508-0004	#508-0847	#508-0005	#508-0848	#508-0338	#508-0863	#508-0864	#508-0865	
W	Steris C-Max	#508-0037	#508-0866	#508-0038	#508-0867	#508-0349	#508-0882	#508-0883	#508-0884	
X	Steris Surgi-Graphic 6000	#508-0885	#508-0886	#508-0887	#508-0888	#508-0905	#508-0906	#508-0907	#508-0908	




	US	SOFTCARE				SOFTCARE PLUS				
		2"	3"	4"	5"	2.25"	3.25"	4.25"	5.25"	
		Metric	5cm	8cm	10cm	13cm	6cm	8cm	11cm	13cm
		PRICE								
A	Berchtold B810, B710	#508-0173	#508-0174	#508-0175	#508-0176	#508-0181	#508-0182	#508-0183	#508-0184	
B	Berchtold Operon D850, D750	#508-0197	#508-0198	#508-0199	#508-0200	#508-0205	#508-0206	#508-0207	#508-0208	
C	Maquet/Siemens 1130	#508-0018	#508-0219	#508-0220	#508-0221	#508-0080	#508-0226	#508-0227	#508-0228	
D	Maquet 1130 Split Leg	#508-0240	#508-0241	#508-0242	#508-0243	#508-0248	#508-0249	#508-0250	#508-0251	
E	Maquet 1130 Orthopedic	#508-0264	#508-0265	#508-0266	#508-0267	#508-0272	#508-0273	#508-0274	#508-0275	
F	Maquet Alphamaxx	#508-0288	#508-0289	#508-0290	#508-0291	#508-0296	#508-0297	#508-0298	#508-0299	
G	Maquet Alphastar 1132	#508-0312	#508-0313	#508-0314	#508-0315	#508-0320	#508-0321	#508-0322	#508-0323	
H	Midmark 7100	#508-0015	#508-0568	#508-0569	#508-0570	#508-0079	#508-0575	#508-0576	#508-0577	
I	Midmark Chic CST	#508-0589	#508-0590	#508-0591	#508-0592	#508-0597	#508-0598	#508-0599	#508-0600	
J	Midmark Montreal	#508-0613	#508-0614	#508-0615	#508-0616	#508-0621	#508-0622	#508-0623	#508-0624	
K	Getinge/Shampaine 1400/1500/1900	#508-0036	#508-1299	#508-1300	#508-1301	#508-0086	#508-1306	#508-1307	#508-1308	
L	Getinge/Shampaine 1501 Radi-Op	#508-0030	#508-0635	#508-0636	#508-0637	#508-0084	#508-0642	#508-0643	#508-0644	
M	Getinge/Shampaine 4900/5100B, 5100E	#508-0033	#508-0654	#508-0655	#508-0656	#508-0085	#508-0661	#508-0662	#508-0663	
N	Skytron 5001, 6001, 6500, 6600	#508-0027	#508-0692	#508-0693	#508-0694	#508-0083	#508-0699	#508-0700	#508-0701	
O	Skytron 3100	#508-0021	#508-0711	#508-0712	#508-0713	#508-0081	#508-0718	#508-0719	#508-0720	
P	Skytron 3500	#508-0024	#508-0730	#508-0731	#508-0732	#508-0082	#508-0737	#508-0738	#508-0739	
Q	Skytron 3501B	#508-0397	#508-0749	#508-0750	#508-0751	#508-0398	#508-0756	#508-0757	#508-0758	
R	Skytron 3600	#508-0770	#508-0771	#508-0772	#508-0773	#508-0778	#508-0779	#508-0780	#508-0781	
S	Skytron 6700, 6701	#508-0392	#508-0792	#508-0793	#508-0794	#508-0393	#508-0799	#508-0800	#508-0801	
T	Steris-Amsco 1080, 2080	#508-0009	#508-0811	#508-0812	#508-0813	#508-0077	#508-0818	#508-0819	#508-0820	
U	Steris-AMSCO 1080/2080 w/Leg Ext.	#508-0012	#508-0830	#508-0831	#508-0832	#508-0078	#508-0837	#508-0838	#508-0839	
V	Steris-Amsco 3080, 3085	#508-0006	#508-0849	#508-0850	#508-0851	#508-0076	#508-0856	#508-0857	#508-0858	
W	Steris C-Max	#508-0039	#508-0868	#508-0869	#508-0870	#508-0087	#508-0875	#508-0876	#508-0877	
X	Steris Surgi-Graphic 6000	#508-0889	#508-0890	#508-0891	#508-0892	#508-0897	#508-0898	#508-0899	#508-0900	



	US	SOFTCARE SILVER				SOFTCARE PLUS SILVER				
		2"	3"	4"	5"	2.25"	3.25"	4.25"	5.25"	
		Metric	5cm	8cm	10cm	13cm	6cm	8cm	11cm	13cm
		PRICE								
A	Berchtold B810, B710	#508-0177	#508-0178	#508-0179	#508-0180	#508-0185	#508-0186	#508-0187	#508-0188	
B	Berchtold Operon D850, D750	#508-0201	#508-0202	#508-0203	#508-0204	#508-0209	#508-0210	#508-0211	#508-0212	
C	Maquet/Siemens 1130	#508-0222	#508-0223	#508-0224	#508-0225	#508-0342	#508-0229	#508-0230	#508-0231	
D	Maquet 1130 Split Leg	#508-0244	#508-0245	#508-0246	#508-0247	#508-0252	#508-0253	#508-0254	#508-0255	
E	Maquet 1130 Orthopedic	#508-0268	#508-0269	#508-0270	#508-0271	#508-0276	#508-0277	#508-0278	#508-0279	
F	Maquet Alphamaxx	#508-0292	#508-0293	#508-0294	#508-0295	#508-0300	#508-0301	#508-0302	#508-0303	
G	Maquet Alphastar 1132	#508-0316	#508-0317	#508-0318	#508-0319	#508-0324	#508-0325	#508-0326	#508-0327	
H	Midmark 7100	#508-0571	#508-0572	#508-0573	#508-0574	#508-0581	#508-0578	#508-0579	#508-0580	
I	Midmark Chic CST	#508-0593	#508-0594	#508-0595	#508-0596	#508-0601	#508-0602	#508-0603	#508-0604	
J	Midmark Montreal	#508-0617	#508-0618	#508-0619	#508-0620	#508-0625	#508-0626	#508-0627	#508-0628	
K	Getinge/Shampaine 1400/1500/1900	#508-1302	#508-1303	#508-1304	#508-1305	#508-1309	#508-1310	#508-1311	#508-1312	
L	Getinge/Shampaine 1501 Radi-Op	#508-0638	#508-0639	#508-0640	#508-0641	#508-0648	#508-0645	#508-0646	#508-0647	
M	Getinge/Shampaine 4900/5100B, 5100E	#508-0657	#508-0658	#508-0659	#508-0660	#508-0667	#508-0664	#508-0665	#508-0666	
N	Skytron 5001, 6001, 6500, 6600	#508-0695	#508-0696	#508-0697	#508-0698	#508-0705	#508-0702	#508-0703	#508-0704	
O	Skytron 3100	#508-0714	#508-0715	#508-0716	#508-0717	#508-0724	#508-0721	#508-0722	#508-0723	
P	Skytron 3500	#508-0733	#508-0734	#508-0735	#508-0736	#508-0743	#508-0740	#508-0741	#508-0742	
Q	Skytron 3501B	#508-0752	#508-0753	#508-0754	#508-0755	#508-0762	#508-0759	#508-0760	#508-0761	
R	Skytron 3600	#508-0774	#508-0775	#508-0776	#508-0777	#508-0782	#508-0783	#508-0784	#508-0785	
S	Skytron 6700, 6701	#508-0795	#508-0796	#508-0797	#508-0798	#508-0805	#508-0802	#508-0803	#508-0804	
T	Steris-Amsco 1080, 2080	#508-0814	#508-0815	#508-0816	#508-0817	#508-0824	#508-0821	#508-0822	#508-0823	
U	Steris-AMSCO 1080/2080 w/Leg Ext.	#508-0833	#508-0834	#508-0835	#508-0836	#508-0843	#508-0840	#508-0841	#508-0842	
V	Steris-Amsco 3080, 3085	#508-0852	#508-0853	#508-0854	#508-0855	#508-0862	#508-0859	#508-0860	#508-0861	
W	Steris C-Max	#508-0871	#508-0872	#508-0873	#508-0874	#508-0881	#508-0878	#508-0879	#508-0880	
X	Steris Surgi-Graphic 6000	#508-0893	#508-0894	#508-0895	#508-0896	#508-0901	#508-0902	#508-0903	#508-0904	



offers a full range of armboard surfaces for any type of armboard in both flat and contoured styles to increase patient arm stability. The Softcare and Softcare Plus Surfaces are designed to reduce pressure and should be part of your pressure management strategy.



	STANDARD		
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)
Allen, Steris-Amsco	#508-0909	#508-0910	#508-0911
Berchtold	#508-0935	#508-0936	#508-0937
Maquet	#508-0961	#508-0962	#508-0963
Midmark	#508-0987	#508-0988	#508-0989
Getinge/Shampaine	#508-1013	#508-1014	#508-1015
Skytron	#508-1039	#508-1040	#508-1041

	SOFTCARE				SOFTCARE PLUS			
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)	5" Pad (13cm)	2.25" Pad (6cm)	3.25" Pad (8cm)	4.25" Pad (11cm)	5.25" Pad (13cm)
Allen, Steris-Amsco	#508-0915	#508-0916	#508-0917	#508-0918	#508-0923	#508-0924	#508-0925	#508-0926
Berchtold	#508-0941	#508-0942	#508-0943	#508-0944	#508-0949	#508-0950	#508-0951	#508-0952
Maquet	#508-0967	#508-0968	#508-0969	#508-0970	#508-0975	#508-0976	#508-0977	#508-0978
Midmark	#508-0993	#508-0994	#508-0995	#508-0996	#508-1001	#508-1002	#508-1003	#508-1004
Getinge/Shampaine	#508-1019	#508-1020	#508-1021	#508-1022	#508-1027	#508-1028	#508-1029	#508-1030
Skytron	#508-1045	#508-1046	#508-1047	#508-1048	#508-1053	#508-1054	#508-1055	#508-1056

	DELUXE			DELUXE GEL			
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)	2.25" Pad (6cm)	3.25" Pad (8cm)	4.25" Pad (11cm)	5.25" Pad (13cm)
Allen, Steris-Amsco	#508-0912	#508-0913	#508-0914	#508-0931	#508-0932	#508-0933	#508-0934
Berchtold	#508-0938	#508-0939	#508-0940	#508-0957	#508-0958	#508-0959	#508-0960
Maquet	#508-0964	#508-0965	#508-0966	#508-0983	#508-0984	#508-0985	#508-0986
Midmark	#508-0990	#508-0991	#508-0992	#508-1009	#508-1010	#508-1011	#508-1012
Getinge/Shampaine	#508-1016	#508-1017	#508-1018	#508-1035	#508-1036	#508-1037	#508-1038
Skytron	#508-1042	#508-1043	#508-1044	#508-1061	#508-1062	#508-1063	#508-1064

	SOFTCARE SILVER				SOFTCARE PLUS SILVER			
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)	5" Pad (13cm)	2.25" Pad (6cm)	3.25" Pad (8cm)	4.25" Pad (11cm)	5.25" Pad (13cm)
Allen, Steris-Amsco	#508-0919	#508-0920	#508-0921	#508-0922	#508-0927	#508-0928	#508-0929	#508-0930
Berchtold	#508-0945	#508-0946	#508-0947	#508-0948	#508-0953	#508-0954	#508-0955	#508-0956
Maquet	#508-0971	#508-0972	#508-0973	#508-0974	#508-0979	#508-0980	#508-0981	#508-0982
Midmark	#508-0997	#508-0998	#508-0999	#508-1000	#508-1005	#508-1006	#508-1007	#508-1008
Getinge/Shampaine	#508-1023	#508-1024	#508-1025	#508-1026	#508-1031	#508-1032	#508-1033	#508-1034
Skytron	#508-1049	#508-1050	#508-1051	#508-1052	#508-1057	#508-1058	#508-1059	#508-1060



ARMBOARD PADS: Contoured Surface

	STANDARD		
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)
Allen, Steris-Amsco	#508-1065	#508-1066	#508-1067
Berchtold	#508-1091	#508-1092	#508-1093
Maquet	#508-1117	#508-1118	#508-1119
Midmark	#508-1143	#508-1144	#508-1145
Getinge/Shampaine	#508-1169	#508-1170	#508-1171
Skytron	#508-1195	#508-1196	#508-1197

	DELUXE		DELUXE GEL				
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)	2.25" Pad (6cm)	3.25" Pad (8cm)	4.25" Pad (11cm)	5.25" Pad (13cm)
Allen, Steris-Amsco	#508-1068	#508-1069	#508-1070	#508-1087	#508-1088	#508-1089	#508-1090
Berchtold	#508-1094	#508-1095	#508-1096	#508-1113	#508-1114	#508-1115	#508-1116
Maquet	#508-1120	#508-1121	#508-1122	#508-1139	#508-1140	#508-1141	#508-1142
Midmark	#508-1146	#508-1147	#508-1148	#508-1165	#508-1166	#508-1167	#508-1168
Getinge/Shampaine	#508-1172	#508-1173	#508-1174	#508-1191	#508-1192	#508-1193	#508-1194
Skytron	#508-1198	#508-1199	#508-1200	#508-1217	#508-1218	#508-1219	#508-1220

ARMBOARD PADS: Contoured Surface



	SOFTCARE				SOFTCARE PLUS			
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)	5" Pad (13cm)	2.25" Pad (6cm)	3.25" Pad (8cm)	4.25" Pad (4.25)	5.25" Pad (13cm)
Allen, Steris-Amsco	#508-1071	#508-1072	#508-1073	#508-1074	#508-1079	#508-1080	#508-1081	#508-1082
Berchtold	#508-1097	#508-1098	#508-1099	#508-1100	#508-1105	#508-1106	#508-1107	#508-1108
Maquet	#508-1123	#508-1124	#508-1125	#508-1126	#508-1131	#508-1132	#508-1133	#508-1134
Midmark	#508-1149	#508-1150	#508-1151	#508-1152	#508-1157	#508-1158	#508-1159	#508-1160
Getinge/Shampaine	#508-1175	#508-1176	#508-1177	#508-1178	#508-1183	#508-1184	#508-1185	#508-1186
Skytron	#508-1201	#508-1202	#508-1203	#508-1204	#508-1209	#508-1210	#508-1211	#508-1212

	SOFTCARE SILVER				SOFTCARE PLUS SILVER			
	2" Pad (5cm)	3" Pad (8cm)	4" Pad (10cm)	5" Pad (13cm)	2.25" Pad (6cm)	3.25" Pad (8cm)	4.25" Pad (11cm)	5.25" Pad (13cm)
Allen, Steris-Amsco	#508-1075	#508-1076	#508-1077	#508-1078	#508-1083	#508-1084	#508-1085	#508-1086
Berchtold	#508-1101	#508-1102	#508-1103	#508-1104	#508-1109	#508-1110	#508-1111	#508-1112
Maquet	#508-1127	#508-1128	#508-1129	#508-1130	#508-1135	#508-1136	#508-1137	#508-1138
Midmark	#508-1153	#508-1154	#508-1155	#508-1156	#508-1161	#508-1162	#508-1163	#508-1164
Getinge/Shampaine	#508-1179	#508-1180	#508-1181	#508-1182	#508-1187	#508-1188	#508-1189	#508-1190
Skytron	#508-1205	#508-1206	#508-1207	#508-1208	#508-1213	#508-1214	#508-1215	#508-1216

SchureMed offers a line of stretcher surfaces with no top edge seams for enhanced infection control. They all come with 4 rounded/contoured corners to fit most stretchers.

	STANDARD					
	74"L (188 cm)	74"L (188 cm)	76"L (193 cm)	76"L (193 cm)	78"L (198 cm)	78"L (198 cm)
	2" (5cm)	3.5" (9cm)	2" (5cm)	3.5" (9cm)	2" (5cm)	3.5" (9cm)
24"	#508-1316	#508-1221	#508-1317	#508-1222	#508-1318	#508-1223
26"	#508-1319	#508-1239	#508-1320	#508-1240	#508-1321	#508-1241
28"	#508-1322	#508-1257	#508-1323	#508-1258	#508-1324	#508-1259
30"	#508-1325	#508-1275	#508-1326	#508-1276	#508-1327	#508-1277

	DELUXE					
	74"L (188 cm)	74"L (188 cm)	76"L (193 cm)	76"L (193 cm)	78"L (98 cm)	78"L (198 cm)
	2" (5cm)	3.5" (9cm)	2" (5cm)	3.5" (9cm)	2" (5cm)	3.5" (9cm)
24"	#508-1328	#508-1224	#508-1329	#508-1225	#508-1330	#508-1226
26"	#508-1331	#508-1242	#508-1332	#508-1243	#508-1333	#508-1244
28"	#508-1334	#508-1260	#508-1335	#508-1261	#508-1336	#508-1262
30"	#508-1337	#508-1278	#508-1338	#508-1279	#508-1339	#508-1280

	SOFTCARE			SOFTCARE PLUS		
	74"L (188 cm)	76"L (193 cm)	78"L (198 cm)	74"L (188 cm)	76"L (193 cm)	78"L (198 cm)
	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)
24"	#508-1227	#508-1228	#508-1229	#508-1233	#508-1234	#508-1235
26"	#508-1245	#508-1246	#508-1247	#508-1251	#508-1252	#508-1253
28"	#508-1263	#508-1264	#508-1265	#508-1269	#508-1270	#508-1271
30"	#508-1281	#508-1282	#508-1283	#508-1287	#508-1288	#508-1289

	SOFTCARE SILVER			SOFTCARE PLUS SILVER		
	74"L (188 cm)	76"L (193 cm)	78"L (198 cm)	74"L (188 cm)	76"L (193 cm)	78"L (198 cm)
	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)	3.5" (9cm)
24"	#508-1230	# 508-1231	#508-1232	#508-1236	#508-1237	#508-1238
26"	#508-1248	#508-1249	#508-1250	#508-1254	#508-1255	#508-1256
28"	#508-1266	#508-1267	#508-1268	#508-1272	#508-1273	#508-1274
30"	#508-1284	#508-1285	#508-1286	#508-1290	#508-1291	#508-1292

Forte de 30 ans d'expériences dans le domaine médical et doué d'un savoir faire unique, notre équipe de professionnels, à la pointe des nouvelles technologies, saura répondre à tous vos besoins: réparations, mise en conformité du matériel, révisions...

Notre réseau de commerciaux formé de manière continue, est à votre service sur l'ensemble du territoire national.



Alphacoméd SARL

ZI Athelia 2
145, impasse du serpolet
13704 La Ciotat Cedex
France

Tel: +33 (0) 04 42 73 10 89
Fax: +33 (0) 04 42 73 10 03
email: alpha@alphacoméd.com
www.alphacoméd.com